



NEW!

HSW MED-CAM ENDOLUX – COMPACT CAMERA SYSTEM WITH INTEGRATED LED LIGHT SOURCE

- Smart workplace concept
- CMOS HD camera
- Endolux LED technology

Compact design and the latest LED technology

The HSW Med-cam Endolux has been developed to provide a compact system for safe diagnosis and surgery. The result is an optimal combination of HD visualization and the latest Endolux LED technology. This leads to bright,

high-resolution images with lifelike colors. Thanks to its compact design and intuitive handling, the system is perfect for the office as well as the operating room.

CMOS HD camera + latest LED technology = HSW Med-cam Endolux



HD camera head with ergonomic design

Endolux LED technology for powerful illumination

All possibilities



HSW Med-cam Endolux	
Article no.	830 002 6454
Illumination	LED
Input voltage range	100–240 VAC, 50/60 Hz
Output	1 x USB for image storage, DVI
HSW camera head	
Article no.	830 001 3118
HSW coupler	
Article no.	830 001 8854

Light guide cable	
Diameter	4.8 mm x 2300 mm
Article no.	830 002 6486

Light guide cable adapters	
Olympus male	Article no. 830 002 6487
ACMI female	Article no. 830 002 9744
Storz/Olympus female	Article no. 830 002 9743
Richard Wolf female	Article no. 830 002 8923