




## ULTRA-SLIM HSW NASOLARYNGOSCOPE

- Ultra-slim with a diameter of 3.4 mm
- Brilliant image quality
- High durability

## Ultra-slim HSW Nasolaryngoscope Premium Endoscopes from Henke-Sass, Wolf

The ultra-slim flexible Nasolaryngoscope 3.4 mm from Henke-Sass, Wolf has been developed in close cooperation with clinical users. The ultra-slim design of the flexible sheath provides an extremely gentle insertion, meaning minimal discomfort for the patients. The excellent optical

quality provides best image for users. The new ultra-slim HSW Nasolaryngoscope impresses with a small diameter, the flexibility of the insertion sheath and an outstanding mechanical stability.



Ultra-slim-line  
for gentle insertion



Precise motion  
for perfect control



Outstanding optical quality  
for brilliant image



High durability  
for cost effectiveness



Ultra-slim HSW Nasolaryngoscope	
Article no.	830 001 9728
Working length	approx. 300 mm
Diameter	3.4 mm
Bending range of distal tip	130° up/down
Direction of view	0°